



Routine cryostat HM 525 UV



Simply
Perfect

Thermo
SCIENTIFIC




Routine cryostat HM 525 UV



The new Microm HM 525 UV takes cryotomy to new levels. With the largest wattage UVC bulb offered in a cryostat, the HM 525 offers the highest level of user security on the market. Stability and performance remain the cornerstone of the HM 525. From the user-friendly controls to the exceptionally stable microtome, the HM 525 UV will make sectioning safer and easier than ever.

- HM 525 UV allows **surface disinfection time from 0-300 minutes**
- 36-Watt UVC bulb offers the largest **wattage** on the market
- Most advanced specimen feed system
- Two independent cutting modes;
fine sectioning 1-20 μm and trim sectioning 10-500 μm



User Security and Easier Sectioning ...

The HM525UVC offers the highest level of user security on the market with a 36 Watt UVC bulb. Fine sectioning 1-20 μm and trim sectioning 10-500 μm .