



## Robot-Stainer HMS 740



**NEW**

### HMS 740

- Automated Documentation of all relevant staining parameters
- Unique Lid-Management saves reagents from oxidization and evaporation
- Drip-Preventing System - less carryover
- Downdraft Ventilation in combination with lids for minimized harmful fumes
- Permanent Independent Agitation for consistent results
- Quality-Management indicates reagent usage

### Technical features:

- 42 overall stations with up to 34 reagent stations
- Touchscreen
- 4 Heat Stations (RT up to 70°C)
- 4 Rinse Stations
- Container: 480ml
- 30 Slides per Basket Carrier
- max. 50 Protocols with 100 steps each

Dimensions: 1160mm x 715mm x 685mm,

Footprint : 1085 x 595mm

Weight:: approx. 98kg



© 2009 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

4481 Campus Drive  
Kalamazoo, MI 49008  
USA  
+1 (800) 522-7270

Robert-Bosch-Strasse 49  
69190 Walldorf  
GERMANY  
+ 49 (0) 6227-8360

[www.thermo.com](http://www.thermo.com)

**Thermo**  
SCIENTIFIC